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
Test

DIAGNOS

FOR IN VITRO DIAGNOSTIC USE ONLY
ONE STEP DIPSTICK TEST

COVID-19 Antigen Test Kit (Colloidal
(Cassette)

For professional in vitro diagnostic use


Please read all the information in the leaflet before performing
MYNORNO USE

The COVID-19 Antigen Test Kit is used for the qualitative detection of protein antigen from SARS-CoV-2 in human throat swabs and nasal swabs. Individuals suspected of COVID-19, only for in vitro diagnostic use. The novel coronavirus belongs to the β genus. COVID-19 is an infectious disease. People are generally susceptible. Currently, the pathogen is the main source of infection. Asymptomatic cases can also be an infectious source. Based on the current epidemiology, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main symptoms include fever, fatigue and dry cough. Nasal congestion, runny nose, muscle pain and diarrhea are found in a few cases.

Antigen is generally detectable in upper respiratory specimens during infection. Rapid diagnosis of SARS-CoV-2 infection will help health-care providers to treat patients and control the disease more efficiently and effectively.

GRAND CLOUTIER

This kit is based on the principle of highly specific antibody-antigen interaction. The test is performed using a colloidal gold labeled immunochromatographic analysis. The test strip contains COVID-19 monoclonal antibody pre-coated in the test area (T) and the COVID-19 monoclonal antibody coated on the label (p). The sample is dropped into the sample well and reacts with the monoclonal antibody which is bound to the pre-coated colloidal gold particles. Then the mixture is chromatographed upwards with capillary action. If the sample is positive, the antibody labeled by colloidal gold particles will first bind to the sample antibody and then to the monoclonal antibody on the test strip.

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A close-up photograph of a person's hands holding a white swab stick. The person is wearing a pink long-sleeved shirt. The swab stick has a white foam tip and a white handle. The person's fingers are wrapped around the handle, and the swab tip is held between their thumb and index finger. The background is a soft, out-of-focus light blue.

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

A woman with dark hair tied back, wearing a white top and a light blue sleeveless vest, is cleaning a light-colored wooden table. She is using a white spray bottle and a white cloth. On the table, there is also a roll of white paper towels. The background shows a modern living room with a white sofa, a bookshelf filled with books, and a large window with sheer white curtains. A green plant is visible on a side table next to the sofa.

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
The box contains : a Sealed pouch, Test Card, Solution Vial, Instructions for use, Swab, Desiccant

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Place the test card on a firm, flat surface. Remove the film from the adhesive strip K

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Disposal of the test kit can be done by adding dishwashing liquid or detergent on the test kit and discarding it in a plastic bag.

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A woman with dark hair in a ponytail, wearing a grey and white sleeveless top, is sitting on a light-colored sofa. She is holding a piece of white trash with Thai text on it. In the background, there is a bookshelf with books, a window with sheer curtains, and a potted plant. A yellow bottle is on a table in the foreground.

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Windows
Go to Settings to activate Windows.

Activate Windows
Go to Settings to activate Windows